
HighRanger 1kb DNA Ladder (300bp - 10,000bp)

Catalogue Number: 11900 100 gel lanes





HighRanger 1kb DNA Ladder (300bp - 10,000bp)Catalogue Number: 11900 (100 gel lanes)

Description:

High quality ready-to-use double-stranded DNA ladder derived by recombinant technology for precise on-gel size and mass determination.

Features:

- Precise band sizes and masses for densitometric quantification.
- Thirteen discrete bands, ranging from 300bp to 10,000bp.
- Higher intensity band at 5,000bp for easy reference.

Contents

1mL of premixed DNA ladder (0.5µg/10µL) in loading buffer (10mM EDTA, 10% glycerol, 0.015% bromophenol blue, and 0.17% SDS)

Instructions:

Mix thoroughly. For best results, load 10μL of DNA ladder per well. For precise mass determination with a densitometer, stain gel after electrophoresis using 0.5μg/ml ethiclium bromide for 30-40 minutes. The accompanying diagram shows the size and mass for each band based on 10μL ladder per well.

Storage:

Stable at room temperature. For longer term storage, -20°C is recommended

This ladder was standardized using $10\mu L$ of DNA per lane on a 0.8 cm thick, 13×15 cm, 1.0% agarose gel run in TAE buffer.

Technical support: techsupport@norgenbiotek.com
Orders: orders@norgenbiotek.com

Toll Free: 1-866-667-4362

Fax: 905-227-1061 www.norgenbiotek.com