

Endotoxin Removal Kit (Mini) - For DNA

Norgen's Endotoxin Removal Kit (Mini) is designed for the rapid removal of endotoxins from up to 25 µg of previously purified DNA. Endotoxins are cell-membrane components of Gram-negative bacteria such as *E. coli*. Endotoxins are released during the lysis step of plasmid purification and significantly reduce transfection efficiencies in endotoxin sensitive cell lines. Therefore, the removal of endotoxins from plasmid preparations is often necessary prior to the use of the DNA in downstream applications. With Norgen's Endotoxin Removal Kit (Mini) endotoxin levels are efficiently reduced to 0.1 EU/ µg DNA or less.



Purification is based on spin column chromatography using Norgen's proprietary resin as the separation matrix. Norgen's resin binds DNA under high salt concentrations and releases the bound DNA under low salt and slightly alkali conditions. The Endotoxin Removal Kit (Mini) can be used to purify plasmids up to 13,000 bp in size, with average recoveries of greater than 90%. The purified DNA is of the highest quality and can be used in a number of downstream applications including sequencing, cloning and transfections.

Kit Specifications			
Maximum DNA Input	25 µg	Size of DNA Purified	Up to 13,000 bp
Maximum DNA Volume Input	100 µL	Average Recovery	> 90%
Final Endotoxin Level	≤ 0.1 EU/µg DNA	Time to Complete 10 Purifications	20 minutes

Endotoxin Removal Kit (Mini) Benefits

Fast and easy processing	Rapid spin-column format allows for the processing of multiple samples in 20 minutes.
Endotoxin-free plasmid DNA	Using Norgen's Endotoxin Removal Kit (Mini), endotoxin levels are reduced to ≤ 0.1 EU / µg of plasmid DNA (Figure 1).
High recovery of input DNA	The average DNA recovery is > 90% of the input amount (Figure 2).
Recovered plasmid DNA is suitable for downstream applications	Purified plasmid DNA is fully compatible with sequencing, cloning and transfections.

Endotoxin Removal Kit (Mini) - For DNA




Figure 1. Efficiently Remove Endotoxins. An input of 11 μg of a plasmid in a 50 μL volume with an endotoxin level of over 8 EU/ μg was subjected to clean-up by Norgen's Endotoxin Removal Kit. One microliter of the input and the 50 μL eluted DNA was diluted 1000-fold or 10 fold, respectively for a 100 μL end-point Limulus Amebocyte Lysate test for endotoxin level. The endotoxin levels of the plasmid DNA sample were reduced from over 8 EU/ μg to less than or equal to 0.1 EU/ μg of plasmid DNA using Norgen's Endotoxin Removal Kit (Mini), with high DNA recovery also (see Figure 2).

Endotoxin Removal Kit (Mini) Contents

1. Binding Solution
2. Wash Solution
3. Elution Buffer
4. Endotoxin Removal Solution
5. Spin Columns
6. Collection Tubes
7. Elution Tubes
8. Product Insert

Shipping Conditions

The Endotoxin Removal Kit (Mini) is shipped at room temperature.




Figure 2. High DNA Recovery. An input of 11 μg of a plasmid in a 50 μL volume with an endotoxin level of over 8 EU/ μg was subjected to clean-up by Norgen's Endotoxin Removal Kit in triplicate. The 50 μL DNA elutions was quantified by spectrophotometry. There is very little loss of DNA associated with using Norgen's Endotoxin Removal Kit (Mini). Recoveries of plasmid DNA are greater than 90% of the input amount when using this kit with efficient endotoxin removal (see Figure 1.)

Customer-Supplied Reagents and Equipment

- Benchtop microcentrifuge
- Microcentrifuge tubes
- 96 - 100% ethanol
- Isopropanol

Storage Conditions

All solutions should be kept tightly sealed and stored at room temperature. All the reagents should remain stable for at least 1 year in their unopened containers.

Cat #	Description	Quantity
22700	Endotoxin Removal Kit (Mini) - For DNA	20 preps