

Bacterial Genomic DNA Isolation Kit

Norgen's Bacterial Genomic DNA Isolation Kit is designed for the rapid preparation of genomic DNA from up to 2×10^9 viable bacterial cells (between 0.5 and 1 mL of culture). Purification is based on spin column chromatography. Norgen's columns bind DNA under high salt concentrations and release the bound DNA under low salt and slightly alkali conditions. The purified DNA is free of RNA and cellular proteinaceous components, and is suitable for many downstream applications.

The Bacterial Genomic DNA Isolation Kit allows for the isolation of genomic DNA from both Gram negative and Gram positive cultures, including *E. coli* and *B. cereus*. Typical purification yields are up to 20 µg of genomic DNA, but will vary depending on the cell density of the bacterial culture and the bacterial species.



Kit Specifications			
Maximum Input	2 x 10 ⁹ bacterial cells	Average Yield*	Up to 20 µg
Column Binding Capacity	25 µg	Time to Complete 10 Purifications	1 hour

* Average yield will vary due to cell density of the bacterial culture, the growth conditions and the bacterial species.

Bacterial Genomic DNA Isolation Kit Benefits

Fast and easy processing	Rapid spin-column format allows for the processing of multiple samples in 1 hour.
High binding capacity of columns	The binding capacity of the columns in the Bacterial Genomic DNA Isolation Kit is 25 µg.
Isolate genomic DNA from all types of bacteria	Genomic DNA can be isolated from both Gram negative and Gram positive bacteria (Figure 1).
Recovered genomic DNA is suitable for downstream applications	Purified genomic DNA is fully compatible with restriction enzyme digestions, sequencing and PCR analysis.
High quality DNA	No degradation of the genomic DNA isolated with the Bacterial Genomic DNA Isolation Kit is observed.

Bacterial Genomic DNA Isolation Kit




Figure 1. Isolation of Genomic DNA from both Gram Positive and Gram Negative Bacteria
The Bacterial Genomic DNA Isolation Kit was used to isolate genomic DNA from the gram-negative bacteria *E. coli* (Lanes 1 and 2), the lysozyme-resistant gram positive bacteria *B. cereus* (Lanes 3 and 4) and the gram positive bacteria *M. luteus* (Lanes 5 and 6). Lane M is Norgen's UltraRanger 1kb DNA Ladder

Bacterial Genomic DNA Isolation Kit Contents

1. Resuspension Buffer
2. Lysis Solution
3. Binding Solution
4. Proteinase K (lyophilized)
5. Wash Solution
6. Elution Buffer
7. Spin Columns inserted into Collection Tubes
8. Elution Tubes
9. Product Insert

Shipping Conditions

The Bacterial Genomic DNA Isolation Kit is shipped at room temperature.




Figure 2. High Yield Purification

The high yield of the Bacterial Genomic DNA Isolation Kit is illustrated by purifying genomic DNA from both a Gram positive and a Gram negative strain, and comparing the yield with two major competitors. With both types of bacteria, Norgen's kit was found to give a higher recovery after both 1 and 2 elutions.

Customer-Supplied Reagents and Equipment

- Benchtop microcentrifuge
- 1.5 mL microcentrifuge tubes
- 55°C water bath or heating block
- 37°C water bath or heating block (gram positive strains only)
- Lysozyme (Gram positive strains only)
- RNase A (optional)

Storage Conditions

All solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for 1 year in their unopened containers. The lyophilized Proteinase K should be stored at -20°C upon arrival and after reconstitution.

Cat #	Description	Quantity
17900	Bacterial Genomic DNA Isolation Kit	50 preps