

NEW Blood/Plasma/Serum RNA, miRNA, DNA & Circulating Nucleic Acid Isolation Kits



Fast, robust, reliable and simple procedures for isolating and purifying high quality, inhibitor-free DNA, RNA, circulating nucleic acids and microRNA from a wide range of sample volumes without phenol-chloroform extraction or alcohol precipitation.

Blood/Plasma/Serum Sample Preparation Kits	Specifications	Catalog Number
Blood Genomic DNA Isolation Micro Kit (50)	Process up to 150 μ L of blood, average yields of 3 - 5 μ g of genomic DNA.	18200
Blood Genomic DNA Isolation Mini Kit (50)	Process from 200 μ L to 500 μ L of blood, average yields of 15 - 30 μ g of genomic DNA.	31000
Blood Genomic DNA Isolation Midi Kit (12)	Process from 500 μ L up to 2 mL of blood, average yields of 30 - 60 μ g of genomic DNA.	31100
Blood Genomic DNA Isolation Maxi Kit (12)	Process up to 2 to 10 mL of blood, average yields of 300 - 600 μ g of genomic DNA.	31200
Plasma/Serum Viral DNA Isolation Kit (50)	Rapid isolation and purification of viral DNA from 0.5 mL plasma or serum samples, compatible with detection products.	29700
Plasma/Serum Circulating DNA Isolation Kit (12)	Rapid concentration of DNA from high sample input into low elution volumes, free of contaminants and inhibitors.	29900
Plasma/Serum Circulating RNA Purification Kit (12)	Rapid concentration of RNA from high sample input into low elution volumes, high efficiency in recovering circulating RNAs both large and small.	30000
Plasma/Serum Circulating RNA Purification 96-Well Kit (Slurry Format) (1 plate)	High throughput concentration of RNA from high sample input to low elution volumes, high efficiency in recovering circulating RNAs of all sizes.	29500
Plasma/Serum Circulating Nucleic Acid Purification Kit (Slurry Format) (25)	Concentrate nucleic acids with high sample input into low elution volumes, high efficiency in recovering circulating DNA and RNA.	27800
Leukocyte RNA Purification Kit (50)	Fractionate leukocytes from whole human blood in minutes, isolate all RNA sizes with an average yield of 1.5 μ g. Compatible with Blood from hamster, mouse, rabbit and rat.	21200

Please contact us for more information regarding these and other innovative sample preparation products from Norgen Biotek—*The Sample Preparation Experts.*

E-mail: info@norgenbiotek.com North American Toll Free: 1-866-667-4362
From Anywhere: 905-227-8848 Fax: 905-227-1061

www.norgenbiotek.com



Norgen Biotek Corp. 3430 Schmon Parkway, Thorold, ON L2V 4Y6