

Complete Line of Urine RNA, miRNA, DNA & Protein Sample Preparation Kits



Fast, robust, reliable and simple procedures for preserving, concentrating, isolating and purifying high quality macromolecules including cell-free total proteins, DNA, RNA and microRNA, from urine.

Product	Specifications	Catalog #
Urine Protein Purification Kits		
Urine Protein Concentration Micro Kit	Concentrate and Desalt 200 - 250 µg of proteins from 0.5 mL to 3 mL of urine within 30 minutes.	17400
Urine Protein Concentration Maxi Kit	Concentrate and Desalt up to 4 mg of proteins from 20 mL of urine within 40 minutes.	21600
Urine DNA/Protein Isolation Kit	Concentrate and Desalt 200 - 250 µg of proteins from 0.5 to 3 mL of urine. Simultaneously isolate exfoliated cell gDNA, apoptotic DNA and urine bacterial DNA within 35 minutes	27400
Urine RNA and microRNA Isolation and Purification		
Urine (Exfoliated Cell) RNA Purification Kit	All RNA species, including microRNA, can be isolated from exfoliated cells found in 1 to 50 mL of urine in 20 minutes.	22500
Urine Bacteria RNA Purification Kit	All RNA species, including microRNA can be isolated from bacterial cells present in 1 to 50 mL of urine in 30 minutes.	23400
Urine microRNA Purification Kit	Isolate and purify small RNA molecules (<200 nt) from 1.5 mL of urine within 20 minutes.	29000
Urine Total RNA Purification Maxi Kit (Slurry Format)	Isolate and purify all RNA species, including microRNA, from 5 to 10 mL of urine within 30 minutes.	29600
Urine DNA Purification Kits		
Urine DNA Isolation Micro Kit	Isolate both high molecule weight DNA and smaller DNA from 2 mL of urine within 35 minutes.	18100
Urine DNA Isolation Mini Kit (Slurry Method)	Isolate both high molecule weight DNA and smaller DNA from 5 mL to 25 mL of urine within 40 minutes.	27000
Urine DNA Isolation Maxi Kit (Slurry Format)	Isolate both high molecule weight DNA and smaller DNA from 20 mL to 25 mL of urine within 40 minutes.	27200
Urine DNA Isolation 96-Well Kit (Slurry Format)	High throughput isolation of both high molecule weight DNA and smaller DNA from 10 mL of urine within 45 minutes	27100
Urine Bacteria DNA Isolation Kit	Isolate genomic DNA from both Gram negative and Gram positive bacteria from 1 to 20 mL of urine within 30 minutes	22400
Urine (Exfoliated Cell) DNA Purification Kit	Isolate genomic DNA from exfoliated cells found in as little as 1 mL of urine within 30 minutes	22300
Urine Collection and Preservation Tubes	Preserve cell-free macromolecules as well as exfoliated cells and microorganisms at ambient temperatures for 1 year	18111

Please contact us for more information regarding these and other innovative sample preparation products from Norgen Biotek—**The Sample Preparation Experts.**

E-mail: info@norgenbiotek.com North American Toll Free: 1-866-667-4362

From Anywhere: 905-227-8848 Fax: 905-227-1061



www.norgenbiotek.com

Norgen Biotek Corp. 3430 Schmon Parkway, Thorold, ON L2V 4Y6